

DYNAMIC FORCE SPECTROSCOPY OF SINGLE MOLECULE LIGAND-RECEPTOR INTERACTIONS AT HIGH LOADING RATES

C.K. Riener, M. Higgins[†], T. Uchihashi[†], S. Yasuda[†], J.E. Sader[‡], R.A. McKendry, S.P. Jarvis[†]
Email: c.riener@ucl.ac.uk

London Centre for Nanotechnology, University College London, WC1E 6JJ, U.K.

[†]Nanoscale Function Group, SFI Nanoscience Laboratory, Trinity College, Dublin 2, Ireland

[‡]Department of Mathematics and Statistics, University of Melbourne, Victoria 3010, Australia

The direct measurement of forces between individual receptor-ligand complexes is an emerging new field. With the atomic force microscope (AFM) it has been possible to probe the forces between a range of biological molecules, including complementary oligo-nucleotides, proteins and carbohydrate interactions. To date, the force required to rupture these single ligand-receptor interactions is commonly detected via static mode deflection, and is typically in the picoNewton range. It has been shown that the measured forces vary with the loading rate, and extrapolation to zero loading rate correlates to the thermal dissociation rate (off-rate) [1].

Here, we introduce the first application of dynamic Frequency Modulated AFM (FM-AFM) to investigate the receptor-ligand interactions at high loading rates. The biotin-avidin interaction is exceptionally strong, $K_D = 10^{-15}$ M, and thus serves as a benchmark for molecular recognition studies. AFM tips were covalently functionalised with biotin using a poly(ethylene glycol) (PEG) crosslinker and avidin was electrostatically immobilised on a mica surface. All measurements were performed in physiological buffer solution [2]. Operating in liquid required the implementation of Magnetically Activated Dynamic mode (MAD mode) [3], whereby a small magnetic particle is attached to the end of a cantilever and an external magnetic field applied via a current carrying coil. Further, by operating the approach in a ‘fly-fishing mode’ direct tip-substrate contact was usually avoided so minimising non-specific interactions.

Frequency shift curves revealed a positive frequency shift in the form of an adhesion peak that corresponded to an unbinding event. The adhesion peak represented the stretching of the PEG linker, subsequent breaking of the bond and further novel rebinding effects. Using the method of Sader and Jarvis [4] to convert frequency shift into force prior to the unbinding instability gave unbinding force values of the order of 300 pN. This is significantly higher than those given in previous studies [2, 5] due to the use of higher loading rates, which in our case correspond to the resonant frequency of the lever (typically of the order of 20 kHz). The significance of these high loading rates will be discussed along with dynamic rebinding effects and the possibility of exploring undiscovered energy barriers in bond-dissociation pathways.

References:

- [1] E. Evans and K. Ritchie (1997). *Biophys. J.* 72:1541–1555
- [2] C.K. Riener et al. (2003). *Anal. Chim. Acta* 479: 59-75
- [3] S. Jarvis et al. (1993). *Rev. Sci. Instrum.* 64: 3515-3520
- [4] J. E. Sader and S. P. Jarvis (2004). *Appl. Phys. Lett.* 84: 1801-1803.
- [5] R. De Paris et al. (2000). *Single Mol.* 4: 285-290